

Anti-Antp antibody (STJA0005709)

STJA0005709

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Antennapedia Homeodomain is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	All

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Antigen Affinity Purified
Dilution Range	WB 1:500-1:1000
Formulation	PBS + 0.02% NaN ₃
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Synthetic peptide corresponding to amino acids from the Antennapedia homeodomain sequence which is 100 % conserved across all known species.
Immunogen Region	
Specificity	Reacts with the Antennapedia homeodomain sequence which is 100 % conserved across all known species.
Immunogen Sequence	