

Anti-MMP9 antibody (STJA0005711)

STJA0005711

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-MMP9 is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Antigen Affinity Purified
Dilution Range	WB 1:500-1:1000
Formulation	PBS + 0.02% NaN ₃
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	4318
Gene Symbol	MMP9
Uniprot ID	MMP9_HUMAN
Immunogen	Synthetic peptide corresponding to amino acids from human MMP-9 protein.
Immunogen Region	
Specificity	Reacts with the 92 kDa human MMP-9 protein.
Immunogen Sequence	