

Anti-F8A1 antibody (C-Term) (STJA0005713)

STJA0005713

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-HAP40 (C-Term) is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Antigen Affinity Purified
Dilution Range	WB 1:500-1:1000
Formulation	PBS + 0.02% NaN ₃
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	474383/474384/8263
Gene Symbol	F8A1.F8A2.F8A3
Uniprot ID	HAP40_HUMAN
Immunogen	Synthetic peptide corresponding to amino acids from the C-terminal region of human HAP40 protein.
Immunogen Region	C-Term
Specificity	Reacts with the C-terminal region of human HAP40 protein.
Immunogen Sequence	