

Anti-KLK4 antibody (C-Term) (STJA0005728)

STJA0005728

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-KLK-L1 (C-Term) is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Antigen Affinity Purified
Dilution Range	WB 1:500-1:1000
Formulation	PBS + 0.02% NaN3
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	9622
Gene Symbol	KLK4
Uniprot ID	KLK4_HUMAN
Immunogen	Synthetic peptide corresponding to amino acids from the C-terminal region of human 39 kDa KLK-L1 (KLK4) protein
Immunogen Region	C-Term
Specificity	Reacts with the C-terminal region of human 39 kDa KLK-L1 (KLK4) protein.
Immunogen Sequence	