

Anti-RETNLA antibody (N-Term) (STJA0005735)

STJA0005735

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-RELM alpha (N-Term) is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Antigen Affinity Purified
Dilution Range	WB 1:500-1:1000
Formulation	PBS + 0.02% NaN ₃
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	57262
Gene Symbol	Retnla
Uniprot ID	RETNA_MOUSE
Immunogen	Synthetic peptide corresponding to amino acids from the N-terminal region of mouse RELM alpha protein
Immunogen Region	N-Term
Specificity	Reacts with the N-terminal region of mouse RELM alpha protein
Immunogen Sequence	