

Anti-BBC3 antibody (STJA0005744)

STJA0005744

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-PUMA is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Antigen Affinity Purified
Dilution Range	WB 1:500-1:1000
Formulation	PBS + 0.02% NaN ₃
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	27113
Gene Symbol	BBC3
Uniprot ID	BBC3_HUMAN
Immunogen	Synthetic peptide corresponding to amino acids from a sequence found in the exon-3-encoded region in the human PUMA-a and PUMA-b forms of the protein.
Immunogen Region	
Specificity	Reacts with a sequence found in the exon-3-encoded region in the human PUMA-a and PUMA-b forms of the protein. Recognizes both forms of the protein
Immunogen Sequence	