

Anti-PERP antibody (First Extracellular Loop) (STJA0005748)

STJA0005748

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-PERP (First Extracellular Loop) is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Antigen Affinity Purified
Dilution Range	WB 1:500-1:1000
Formulation	PBS + 0.02% NaN ₃
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	64065
Gene Symbol	PERP
Uniprot ID	PERP_HUMAN
Immunogen	Synthetic peptide corresponding to amino acids from the first extracellular loop region of the human PERP protein.
Immunogen Region	First Extracellular Loop
Specificity	Reacts with the first extracellular loop region of the human PERP protein.
Immunogen Sequence	