

Anti-COL25A1 antibody (NC3 region) (STJA0005756)

STJA0005756

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-CLAC-P (NC3 region) is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Antigen Affinity Purified
Dilution Range	WB 1:500-1:1000
Formulation	PBS + 0.02% NaN ₃
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	84570
Gene Symbol	COL25A1
Uniprot ID	COPA1_HUMAN
Immunogen	Synthetic peptide corresponding to amino acids from the NC23 region of human CLAC-P
Immunogen Region	NC3 region
Specificity	Reacts with the NC3 region of human and mouse CLAC-P.
Immunogen Sequence	