

**Anti-RASSF5 antibody (STJA0005761)**

STJA0005761

**GENERAL INFORMATION**

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Rabbit polyclonal antibody anti-RASSF1 is suitable for use in Western Blot research applications.
<b>Applications</b>	WB
<b>Host/Source</b>	Rabbit
<b>Reactivity</b>	Human/Mouse/Rat

**PRODUCT PROPERTIES**

<b>Clonality</b>	Polyclonal
<b>Clone ID</b>	
<b>Concentration</b>	
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	Antigen Affinity Purified
<b>Dilution Range</b>	WB 1:500-1:1000
<b>Formulation</b>	PBS + 0.02% NaN <sub>3</sub>
<b>Isotype</b>	IgG
<b>Storage Instruction</b>	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

**TARGET INFORMATION**

<b>Gene ID</b>	<a href="#">83593</a>
<b>Gene Symbol</b>	<a href="#">RASSF5</a>
<b>Uniprot ID</b>	<a href="#">RASF5_HUMAN</a>
<b>Immunogen</b>	Synthetic peptide corresponding to amino acids from the RA domain of human NORE1
<b>Immunogen Region</b>	
<b>Specificity</b>	Reacts with residues located in the RA domain of human, mouse and rat NORE1 protein.
<b>Immunogen Sequence</b>	