

Anti-MINDY3 antibody (N-Term) (STJA0005762)

STJA0005762

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-CARP (N-Term) is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Antigen Affinity Purified
Dilution Range	WB 1:500-1:1000
Formulation	PBS + 0.02% NaN ₃
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	80013
Gene Symbol	MINDY3
Uniprot ID	MINDY3_HUMAN
Immunogen	Synthetic peptide corresponding to amino acids from the N-terminal region of human CARP protein
Immunogen Region	N-Term
Specificity	Reacts with the N-terminal region of human CARP protein.
Immunogen Sequence	