

Anti-PRG3 antibody (STJA0005763)

STJA0005763

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-PRG3 is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Antigen Affinity Purified
Dilution Range	WB 1:500-1:1000
Formulation	PBS + 0.02% NaN ₃
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	10394
Gene Symbol	PRG3
Uniprot ID	PRG3_HUMAN
Immunogen	
Immunogen Region	
Specificity	Reacts with the 41 kDa human PRG3 protein. This sequence is 87% identical in mouse.
Immunogen Sequence	