

**Anti-NSMCE2 antibody (C-Term) (STJA0005768)**

STJA0005768

**GENERAL INFORMATION**

|                          |                                                                                                          |
|--------------------------|----------------------------------------------------------------------------------------------------------|
| <b>Product Type</b>      | Primary antibodies                                                                                       |
| <b>Short Description</b> | Rabbit polyclonal antibody anti-NSE2 (C-Term) is suitable for use in Western Blot research applications. |
| <b>Applications</b>      | WB                                                                                                       |
| <b>Host/Source</b>       | Rabbit                                                                                                   |
| <b>Reactivity</b>        | Human                                                                                                    |

**PRODUCT PROPERTIES**

|                            |                                                                                                |
|----------------------------|------------------------------------------------------------------------------------------------|
| <b>Clonality</b>           | Polyclonal                                                                                     |
| <b>Clone ID</b>            |                                                                                                |
| <b>Concentration</b>       |                                                                                                |
| <b>Conjugation</b>         | Unconjugated                                                                                   |
| <b>Purification</b>        | Antigen Affinity Purified                                                                      |
| <b>Dilution Range</b>      | WB 1:500-1:1000                                                                                |
| <b>Formulation</b>         | PBS + 0.02% NaN <sub>3</sub>                                                                   |
| <b>Isotype</b>             | IgG                                                                                            |
| <b>Storage Instruction</b> | Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles. |

**TARGET INFORMATION**

|                           |                                                                                          |
|---------------------------|------------------------------------------------------------------------------------------|
| <b>Gene ID</b>            | <a href="#">286053</a>                                                                   |
| <b>Gene Symbol</b>        | <a href="#">NSMCE2</a>                                                                   |
| <b>Uniprot ID</b>         | <a href="#">NSE2_HUMAN</a>                                                               |
| <b>Immunogen</b>          | Synthetic peptide corresponding to amino acids from the C-terminal region of human NSE2. |
| <b>Immunogen Region</b>   | C-Term                                                                                   |
| <b>Specificity</b>        | Reacts with the C-terminal region of human Nse2                                          |
| <b>Immunogen Sequence</b> |                                                                                          |