

Anti-FERMT1 antibody (C-Term) (STJA0005772)

STJA0005772

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-URP1 (C-Term) is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Bovine/Canine/Human/Mouse/Non-Human Primates/Rat/Xenopus

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Antigen Affinity Purified
Dilution Range	WB 1:500-1:1000
Formulation	PBS + 0.02% NaN ₃
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	55612
Gene Symbol	FERMT1
Uniprot ID	FERM1_HUMAN
Immunogen	Synthetic peptide corresponding to amino acids from the C-terminal region of human URP1.
Immunogen Region	C-Term
Specificity	Reacts with the C-terminal region of the 77 kDa human URP1 protein. Sequence is 100% conserved between almost all species (mouse, rat, canine, bovine, chimpanzee, Xenopus).
Immunogen Sequence	