

Anti-XRCC5 antibody (Internal Region) (STJA0005773)

STJA0005773

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Ku80 (Internal Region) is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Hamster/Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Antigen Affinity Purified
Dilution Range	WB 1:500-1:1000
Formulation	PBS + 0.02% NaN ₃
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	7520
Gene Symbol	XRCC5
Uniprot ID	XRCC5_HUMAN
Immunogen	Synthetic peptide corresponding to amino acids from an internal sequence of the human Ku80 protein
Immunogen Region	Internal Region
Specificity	Reacts with the 80 kDa Ku80 protein. Sequence is 100% conserved between human, mouse, rat and hamster.
Immunogen Sequence	