

**Anti-XRCC5 antibody (Internal Region) (STJA0005774)**

STJA0005774

**GENERAL INFORMATION**

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Rabbit polyclonal antibody anti-Ku80 (Internal Region) is suitable for use in Western Blot research applications.
<b>Applications</b>	WB
<b>Host/Source</b>	Rabbit
<b>Reactivity</b>	Hamster/Human/Mouse/Rat

**PRODUCT PROPERTIES**

<b>Clonality</b>	Polyclonal
<b>Clone ID</b>	
<b>Concentration</b>	
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	Antigen Affinity Purified
<b>Dilution Range</b>	WB 1:500-1:1000
<b>Formulation</b>	PBS + 0.02% NaN <sub>3</sub>
<b>Isotype</b>	IgG
<b>Storage Instruction</b>	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

**TARGET INFORMATION**

<b>Gene ID</b>	<a href="#">7520</a>
<b>Gene Symbol</b>	<a href="#">XRCC5</a>
<b>Uniprot ID</b>	<a href="#">XRCC5_HUMAN</a>
<b>Immunogen</b>	Synthetic peptide corresponding to amino acids from an internal sequence of the human Ku80 protein
<b>Immunogen Region</b>	Internal Region
<b>Specificity</b>	Reacts with the 80 kDa Ku80 protein. Sequence is 100% conserved between human, mouse, rat and hamster.
<b>Immunogen Sequence</b>	