

Anti-ZNRF2 antibody (Internal Region) (STJA0005780)

STJA0005780

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-ZNRF2 (Internal Region) is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Antigen Affinity Purified
Dilution Range	WB 1:500-1:1000
Formulation	PBS + 0.02% NaN ₃
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	223082
Gene Symbol	ZNRF2
Uniprot ID	ZNRF2_HUMAN
Immunogen	Synthetic peptide corresponding to amino acids from an internal region of human DPP10 protein
Immunogen Region	Internal Region
Specificity	Reacts with an internal region of the 24 kDa human ZNRF2 protein. 87% sequence identity to rat and 81% sequence identity to mouse.
Immunogen Sequence	