

Anti-LY6K antibody (Internal Region) (STJA0005781)

STJA0005781

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Ly-6k (Internal Region) is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Antigen Affinity Purified
Dilution Range	WB 1:500-1:1000
Formulation	PBS + 0.02% NaN ₃
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	54742
Gene Symbol	LY6K
Uniprot ID	LY6K_HUMAN
Immunogen	Synthetic peptide corresponding to amino acids from an internal sequence of the human Ly-6k protein
Immunogen Region	Internal Region
Specificity	Reacts with the 14 kDa human Ly-6K protein
Immunogen Sequence	