

Anti-BCAS1 antibody (N-Term) (STJA0005783)

STJA0005783

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-NABC1 (N-Term) is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Antigen Affinity Purified
Dilution Range	WB 1:500-1:1000
Formulation	PBS + 0.02% NaN ₃
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	8537
Gene Symbol	BCAS1
Uniprot ID	BCAS1_HUMAN
Immunogen	Synthetic peptide corresponding to amino acids from the N-terminal region of human NABC1 protein
Immunogen Region	N-Term
Specificity	Reacts with the C-terminal region of human NABC1. NABC1 is a 60 kDa protein; however it forms dimers, which are easily detected in MCF-7 cells.
Immunogen Sequence	