

Anti-MSI1 antibody (N-Term) (STJA0005789)

STJA0005789

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Musashi 1 (N-Term) is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Antigen Affinity Purified
Dilution Range	WB 1:500-1:1000
Formulation	PBS + 0.02% NaN ₃
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	4440
Gene Symbol	MSI1
Uniprot ID	MSI1H_HUMAN
Immunogen	Synthetic peptide corresponding to amino acids from the N-terminal region of the human, mouse and rat Musashi 1 protein
Immunogen Region	N-Term
Specificity	Reacts with the N-terminal region of the human, mouse and rat Musashi 1 protein
Immunogen Sequence	