

Anti-CD9 antibody (N-Term) (STJA0005794)

STJA0005794

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-CD9 (N-Term) is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Antigen Affinity Purified
Dilution Range	WB 1:500-1:1000
Formulation	PBS + 0.02% NaN ₃
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	928
Gene Symbol	CD9
Uniprot ID	CD9_HUMAN
Immunogen	Synthetic peptide corresponding to amino acids from the N-terminal region of human CD9 protein.
Immunogen Region	N-Term
Specificity	Reacts with the N-terminal region of the 25 kDa human, mouse and rat CD9 protein.
Immunogen Sequence	