

Anti-TUBB3 antibody (C-Term) (STJA0005804)

STJA0005804

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Beta Tubulin III (C-Term) is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Antigen Affinity Purified
Dilution Range	WB 1:500-1:1000
Formulation	PBS + 0.02% NaN ₃
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	10381
Gene Symbol	TUBB3
Uniprot ID	TBB3_HUMAN
Immunogen	Synthetic peptide corresponding to amino acids from the C-terminal region human beta-tubulin-III protein
Immunogen Region	C-Term
Specificity	Reacts with C-terminal region of the 50 kDa human beta-tubulin-III protein. 94% sequence identity in mouse and rat.
Immunogen Sequence	