

Anti-NANOG antibody (Internal Region) (STJA0005813)

STJA0005813

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Nanog (Internal Region) is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Antigen Affinity Purified
Dilution Range	WB 1:500-1:1000
Formulation	PBS + 0.02% NaN ₃
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	79923
Gene Symbol	NANOG
Uniprot ID	NANOG_HUMAN
Immunogen	Synthetic peptide corresponding to an internal sequence of the human Nanog protein.
Immunogen Region	Internal Region
Specificity	Reacts with an internal sequence of the human nanog protein
Immunogen Sequence	