

Anti-Cleaved-F11-Arg387, Heavy Chain antibody (387-388) (STJA0005836)

STJA0005836

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Cleaved-Coagulation factor XI-Arg387, Heavy Chain (387-388) is suitable for use in Western Blot and ELISA research applications.
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution Range	WB 1:1000-2000 ELISA 1:5000-20000
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	2160
Gene Symbol	F11
Uniprot ID	FA11_HUMAN
Immunogen	Synthesized peptide derived from human FA11 (heavy chain, Cleaved-Arg387)
Immunogen Region	387-388
Specificity	This antibody detects endogenous levels of Human, Mouse FA11 (heavy chain, Cleaved-Arg387, protein was cleaved amino acid sequence between 387-388)
Immunogen Sequence	