

Anti-Cleaved-COL18A1-His1572 antibody (1571-1572) (STJA0005845)

STJA0005845

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Cleaved-Collagen alpha-1XVIII chain-His1572 (1571-1572) is suitable for use in Western Blot and ELISA research applications.
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution Range	WB 1:1000-2000 ELISA 1:5000-20000
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	80781
Gene Symbol	COL18A1
Uniprot ID	COIA1_HUMAN
Immunogen	Synthesized peptide derived from human Collagen XVIII Alpha 1 (Cleaved-His1572)
Immunogen Region	1571-1572
Specificity	This antibody detects endogenous levels of Human, Mouse Collagen XVIII Alpha 1 (Cleaved-His1572, protein was cleaved amino acid sequence between 1571-1572)
Immunogen Sequence	