

**Anti-Cleaved-COL4A2-Ser1485 antibody (1485-1486) (STJA0005846)**

STJA0005846

**GENERAL INFORMATION**

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Rabbit polyclonal antibody anti-Cleaved-Collagen alpha-2IV chain-Ser1485 (1485-1486) is suitable for use in Western Blot and ELISA research applications.
<b>Applications</b>	WB/ELISA
<b>Host/Source</b>	Rabbit
<b>Reactivity</b>	Human/Mouse

**PRODUCT PROPERTIES**

<b>Clonality</b>	Polyclonal
<b>Clone ID</b>	
<b>Concentration</b>	1 mg/mL
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
<b>Dilution Range</b>	WB 1:1000-2000 ELISA 1:5000-20000
<b>Formulation</b>	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
<b>Isotype</b>	IgG
<b>Storage</b>	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
<b>Instruction</b>	

**TARGET INFORMATION**

<b>Gene ID</b>	<a href="#">1284</a>
<b>Gene Symbol</b>	<a href="#">COL4A2</a>
<b>Uniprot ID</b>	<a href="#">CO4A2_HUMAN</a>
<b>Immunogen</b>	Synthesized peptide derived from human Collagen IV Alpha 2 (Cleaved-Ser1485)
<b>Immunogen Region</b>	1485-1486
<b>Specificity</b>	This antibody detects endogenous levels of Human, Mouse Collagen IV Alpha 2 (Cleaved-Ser1485, protein was cleaved amino acid sequence between 1485-1486 )
<b>Immunogen Sequence</b>	