

Anti-Cleaved-TMPRSS15-Ile785, Light Chain antibody (785-786) (STJA0005856)

STJA0005856

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Cleaved-Enteropeptidase-Ile785, Light Chain (785-786) is suitable for use in Western Blot and ELISA research applications.
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Rat/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution Range	WB 1:1000-2000 ELISA 1:5000-20000
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	5651
Gene Symbol	TMPRSS15
Uniprot ID	ENTK_HUMAN
Immunogen	Synthesized peptide derived from human ENTK (light chain, Cleaved-Ile785)
Immunogen Region	785-786
Specificity	This antibody detects endogenous levels of Human ENTK (light chain, Cleaved-Ile785, protein was cleaved amino acid sequence between 785-786)
Immunogen Sequence	