

Anti-Cleaved-EFNA1-Ser182 antibody (182-183) (STJA0005858)

STJA0005858

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Cleaved-Ephrin-A1-Ser182 (182-183) is suitable for use in Western Blot and ELISA research applications.
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution Range	WB 1:1000-2000 ELISA 1:5000-20000
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	1942
Gene Symbol	EFNA1
Uniprot ID	EFNA1_HUMAN
Immunogen	Synthesized peptide derived from human EFNA1 (Cleaved-Ser182)
Immunogen Region	182-183
Specificity	This antibody detects endogenous levels of Human, Mouse, Rat EFNA1 (Cleaved-Ser182, protein was cleaved amino acid sequence between 182-183)
Immunogen Sequence	