

## Anti-Cleaved-GGT5-Gly387, Heavy Chain antibody (387-388) (STJA0005867)

STJA0005867

### GENERAL INFORMATION

|                          |                                                                                                                                                                              |
|--------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Type</b>      | Primary antibodies                                                                                                                                                           |
| <b>Short Description</b> | Rabbit polyclonal antibody anti-Cleaved-Glutathione hydrolase 5 proenzyme-Gly387, Heavy Chain (387-388) is suitable for use in Western Blot and ELISA research applications. |
| <b>Applications</b>      | WB/ELISA                                                                                                                                                                     |
| <b>Host/Source</b>       | Rabbit                                                                                                                                                                       |
| <b>Reactivity</b>        | Human/Mouse/Rat                                                                                                                                                              |

### PRODUCT PROPERTIES

|                       |                                                                                                           |
|-----------------------|-----------------------------------------------------------------------------------------------------------|
| <b>Clonality</b>      | Polyclonal                                                                                                |
| <b>Clone ID</b>       |                                                                                                           |
| <b>Concentration</b>  | 1 mg/mL                                                                                                   |
| <b>Conjugation</b>    | Unconjugated                                                                                              |
| <b>Purification</b>   | The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen. |
| <b>Dilution Range</b> | WB 1:1000-2000<br>ELISA 1:5000-20000                                                                      |
| <b>Formulation</b>    | Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.                                   |
| <b>Isotype</b>        | IgG                                                                                                       |
| <b>Storage</b>        | Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.            |
| <b>Instruction</b>    |                                                                                                           |

### TARGET INFORMATION

|                           |                                                                                                                                                           |
|---------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Gene ID</b>            | <a href="#">2687</a>                                                                                                                                      |
| <b>Gene Symbol</b>        | <a href="#">GGT5</a>                                                                                                                                      |
| <b>Uniprot ID</b>         | <a href="#">GGT5_HUMAN</a>                                                                                                                                |
| <b>Immunogen</b>          | Synthesized peptide derived from human GGT5 (heavy chain, Cleaved-Gly387)                                                                                 |
| <b>Immunogen Region</b>   | 387-388                                                                                                                                                   |
| <b>Specificity</b>        | This antibody detects endogenous levels of Human, Mouse, Rat GGT5 (heavy chain, Cleaved-Gly387, protein was cleaved amino acid sequence between 387-388 ) |
| <b>Immunogen Sequence</b> |                                                                                                                                                           |