

Anti-Cleaved-GGT1-Gly380, Heavy Chain antibody (380-381) (STJA0005868)

STJA0005868

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Cleaved-Glutathione hydrolase 1 proenzyme-Gly380, Heavy Chain (380-381) is suitable for use in Western Blot and ELISA research applications.
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Rat/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution Range	WB 1:1000-2000 ELISA 1:5000-20000
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	2678
Gene Symbol	GGT1
Uniprot ID	GGT1_HUMAN
Immunogen	Synthesized peptide derived from human GGT1 (heavy chain, Cleaved-Gly380)
Immunogen Region	380-381
Specificity	This antibody detects endogenous levels of Human GGT1 (heavy chain, Cleaved-Gly380, protein was cleaved amino acid sequence between 380-381)
Immunogen Sequence	