

Anti-Cleaved-KLK1/2/3-Ile25 antibody (24-25) (STJA0005883)

STJA0005883

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Cleaved-KLK1/2/3-Ile25 (24-25) is suitable for use in Western Blot and ELISA research applications.
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Rat/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution Range	WB 1:1000-2000 ELISA 1:5000-20000
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	3816 3817 354 KLK1 KLK2 KLK1_HUMAN KLK2_HUMAN KLK3_HUMAN
Immunogen	Synthesized peptide derived from human KLK1/2/3 (Cleaved-Ile25)
Immunogen Region	24-25
Specificity	This antibody detects endogenous levels of Human KLK1/2/3 (Cleaved-Ile25, protein was cleaved amino acid sequence between 24-25)
Immunogen Sequence	