

Anti-Cleaved-MCFD2-Glu27 antibody (26-27) (STJA0005899)

STJA0005899

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Cleaved-Multiple coagulation factor deficiency protein 2-Glu27 (26-27) is suitable for use in Western Blot and Immunohistochemistry research applications.
Applications	WB/IHC
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution Range	WB 1:500-2000 IHC-p 1:50-300
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	90411
Gene Symbol	MCFD2
Uniprot ID	MCFD2_HUMAN
Immunogen	Synthesized peptide derived from human MCFD2 (Cleaved-Glu27)
Immunogen Region	26-27
Specificity	This antibody detects endogenous levels of Human, Mouse, Rat MCFD2 (Cleaved-Glu27, protein was cleaved amino acid sequence between 26-27)
Immunogen Sequence	



Immunohistochemical analysis of paraffin-embedded human liver cancer. 1. Antibody was diluted at 1:200 (4A°C overnight). 2. Tris-EDTA, pH9.0 was used for antigen retrieval. 3. Secondary antibody was diluted at 1:200 (room temperature, 45min).