

Anti-Acetyl-ACSS2-Lys418 antibody (STJA0005918)

STJA0005918

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Acetyl-Acetyl-coenzyme A synthetase, cytoplasmic-Lys418 is suitable for use in Western Blot and ELISA research applications.
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution Range	WB 1:1000-2000 ELISA 1:5000-20000
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	55902
Gene Symbol	ACSS2
Uniprot ID	ACSA_HUMAN
Immunogen	Synthesized peptide derived from human ACSS2 (Acetyl Lys418)
Immunogen Region	
Specificity	This antibody detects endogenous levels of Human, Mouse ACSS2 (Acetyl Lys418)
Immunogen Sequence	