

Anti-Acetyl-H2AC11; H2AC13; H2AC15; H2AC16; H2AC17-Lys6 antibody (STJA0005937)
STJA0005937**GENERAL INFORMATION**

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Acetyl-Histone H2A type 1-Lys6 is suitable for use in Western Blot and ELISA research applications.
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution Range	WB 1:1000-2000 ELISA 1:5000-20000
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	8329/8330/8332/8336/8969
Gene Symbol	H2AC11.H2AC13.H2AC15.H2AC16.H2AC17
Uniprot ID	H2A1_HUMAN
Immunogen	Synthesized peptide derived from human Histone H2A (Acetyl Lys6)
Immunogen Region	
Specificity	This antibody detects endogenous levels of Human, Mouse, Rat Histone H2A (Acetyl Lys6)
Immunogen Sequence	