

Anti-Acetyl-EEF1A1-Lys146 antibody (STJA0005952)

STJA0005952

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Acetyl-EF-1 Alpha 1/2-Lys146 is suitable for use in Western Blot and ELISA research applications.
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution Range	WB 1:1000-2000 ELISA 1:5000-20000
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	1917 NA 1915 1917 NA
Gene Symbol	EEF1A2 EEF1A1P5 EF1A2_HUMAN EF1A3_HUMAN EF1A1_HUMAN
Immunogen	Synthesized peptide derived from human EF-1 Alpha 1/2 (Acetyl Lys146)
Immunogen Region	
Specificity	This antibody detects endogenous levels of Human, Mouse, Rat EF-1 Alpha 1/2 (Acetyl Lys146)
Immunogen Sequence	