

Anti-Pan-Acetyl-KRT2-Lys185 antibody (STJA0005971)
STJA0005971

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Pan-Acetyl-Keratin-Lys185 is suitable for use in Western Blot and ELISA research applications.
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution Range	WB 1:1000-2000 ELISA 1:5000-20000
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	9119 3856 3849 3854 51350 3853 3890 3850 286887 121391 319101 112802 338785 140807 3852 3855 9119 3856 3849 3854 51350 3853 3890 3850 286887 121391 319101 112802 338785 140807 3852 3855 9119 3856 3849 3854 51350 3853 3890 3850 286887 121391 319101 11
Gene Symbol	KRT75 KRT8 KRT2 KRT6B KRT76 KRT6A KRT84 KRT3 KRT6C KRT74 KRT73 KRT71 KRT79 KRT72 KRT5 KRT7 KRT75 KRT8 KRT2 KRT6B KRT76 KRT6A KRT84 KRT3 KRT6C KRT74 KRT73 KRT71 KRT79 KRT72 KRT5 KRT7 KRT75 KRT8 KRT2 KRT6B KRT76 KRT6A KRT84 KRT3 KRT6C KRT74 KRT73 KRT71
Uniprot ID	K2C75_HUMAN K2C8_HUMAN K22E_HUMAN K2C6B_HUMAN K22O_HUMAN K2C6A_HUMAN KRT84_HUMAN K2C3_HUMAN K2C6C_HUMAN K2C74_HUMAN K2C73_HUMAN K2C71_HUMAN K2C79_HUMAN K2C72_HUMAN K2C5_HUMAN K2C7_HUMAN K2C75_HUMAN K2C8_HUMAN K22E_HUMAN K2C6B_HUMAN K22O_HUMAN K2C6A_H
Immunogen	Synthesized peptide derived from human Keratin-pan (Acetyl Lys185)
Immunogen Region	
Specificity	This antibody detects endogenous levels of Human, Mouse, Rat Keratin-pan (Acetyl Lys185)
Immunogen Sequence	