

Anti-Acetyl-POLB-Lys72 antibody (STJA0005987)

STJA0005987

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Acetyl-DNA polymerase beta-Lys72 is suitable for use in Western Blot and ELISA research applications.
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution Range	WB 1:1000-2000 ELISA 1:5000-20000
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	5423
Gene Symbol	POLB
Uniprot ID	DPOLB_HUMAN
Immunogen	Synthesized peptide derived from human DNA pol Beta (Acetyl Lys72)
Immunogen Region	
Specificity	This antibody detects endogenous levels of Human, Mouse, Rat DNA pol Beta (Acetyl Lys72)
Immunogen Sequence	