

## Anti-Acetyl-TUBA1A-Lys163 antibody (STJA0005998)

STJA0005998

### GENERAL INFORMATION

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Rabbit polyclonal antibody anti-Acetyl-Tubulin alpha-Lys163 is suitable for use in Western Blot and ELISA research applications.
<b>Applications</b>	WB/ELISA
<b>Host/Source</b>	Rabbit
<b>Reactivity</b>	Human/Mouse/Rat

### PRODUCT PROPERTIES

<b>Clonality</b>	Polyclonal
<b>Clone ID</b>	
<b>Concentration</b>	1 mg/mL
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
<b>Dilution Range</b>	WB 1:1000-2000 ELISA 1:5000-20000
<b>Formulation</b>	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
<b>Isotype</b>	IgG
<b>Storage Instruction</b>	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

### TARGET INFORMATION

<b>Gene ID</b>	<a href="#">10376</a> <a href="#">7277</a> <a href="#">7846</a> <a href="#">TUBA1B</a> <a href="#">TUBA4A</a> <a href="#">TBA1B_HUMAN</a> <a href="#">TBA4A_HUMAN</a> <a href="#">TBA1A_HUMAN</a>
<b>Immunogen</b>	Synthesized peptide derived from human Tubulin Alpha (Acetyl Lys163)
<b>Immunogen Region</b>	
<b>Specificity</b>	This antibody detects endogenous levels of Human, Mouse, Rat Tubulin Alpha (Acetyl Lys163)
<b>Immunogen Sequence</b>	