

Anti-Acetyl-TUBB2A-Lys379 antibody (STJA0005999)

STJA0005999

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Acetyl-Tubulin beta-Lys379 is suitable for use in Western Blot and ELISA research applications.
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution Range	WB 1:1000-2000 ELISA 1:5000-20000
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	347733 10381 84617
Gene Symbol	TUBB2B TUBB3 TBB2B_HUMAN TBB3_HUMAN TBB6_HUMAN
Immunogen	Synthesized peptide derived from human Tubulin Beta (Acetyl Lys379)
Immunogen Region	
Specificity	This antibody detects endogenous levels of Human, Mouse, Rat Tubulin Beta (Acetyl Lys379)
Immunogen Sequence	