

Anti-Acetyl-TUBA1A-Lys401 antibody (STJA0006000)

STJA0006000

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Acetyl-Tubulin alpha-Lys401 is suitable for use in Western Blot and ELISA research applications.
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution Range	WB 1:1000-2000 ELISA 1:5000-20000
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	7846 84790 10376 TUBA1A TUBA1C TBA1A_HUMAN TBA1C_HUMAN TBA1B_HUMAN
Immunogen	Synthesized peptide derived from human Tubulin Alpha (Acetyl Lys401)
Immunogen Region	
Specificity	This antibody detects endogenous levels of Human, Mouse, Rat Tubulin Alpha (Acetyl Lys401)
Immunogen Sequence	