

Anti-Pan-Acetyl-ACTA2-Lys63/61/62 antibody (STJA0006005)

STJA0006005

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Pan-Acetyl-Actin-pan-Lys63/61/62 is suitable for use in Western Blot and ELISA research applications.
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution Range	WB 1:1000-2000 ELISA 1:5000-20000
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	58 71 60 ACTA1 ACTG1 ACTS_HUMAN ACTG_HUMAN ACTB_HUMAN
Immunogen	Synthesized peptide derived from human Actin-pan (Acetyl Lys63/61/62)
Immunogen Region	
Specificity	This antibody detects endogenous levels of Human, Mouse, Rat Actin-pan (Acetyl Lys63/61/62)
Immunogen Sequence	