

Anti-PPIB antibody [2B10] {AbFluor™ 405} (STJA0006027)

STJA0006027

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Peptidyl-prolyl cis-trans isomerase B is suitable for use in Immunofluorescence, Western Blot and Immunohistochemistry research applications.
Applications	IF/WB/IHC
Host/Source	Mouse
Reactivity	Human/Rat/Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	2B10
Concentration	1 mg/mL
Conjugation	AbFluor™ 405
Purification	The antibody was purified using affinity-chromatography using specific immunogen.
Dilution Range	Optimal working dilutions should be determined experimentally by the investigator Suggested starting dilutions are as follows: IHC 1:50-300 IF 1:50-200
Formulation	Liquid in PBS pH7.4, containing 0.02% Sodium Azide and 50% Glycerol.
Isotype	IgG
Storage	Stable for one year at -15°C to -25°C from date of shipment. For maximum recovery of product, centrifuge the original vial after thawing and prior to removing the cap. Aliquot to avoid repeated freezing
Instruction	

TARGET INFORMATION

Gene ID	5479
Gene Symbol	PPIB
Uniprot ID	PPIB_HUMAN
Immunogen	
Immunogen Region	
Specificity	Cyclophilin B Monoclonal Antibody (2B10) AbFluor™ 405 Conjugated specially designed for your Immunofluorescence analysis.
Immunogen Sequence	