

Anti-Flag-Tag antibody [2C5] {AbFluor™ 405} (STJA0006039)

STJA0006039

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Flag-Tag-Tag is suitable for use in Western Blot, ELISA, Immunoprecipitation and Immunofluorescence research applications.
Applications	WB/ELISA/IP/IF
Host/Source	Mouse
Reactivity	Species independent

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	2C5
Concentration	1 mg/mL
Conjugation	AbFluor™ 405
Purification	The antibody was purified using affinity-chromatography using specific immunogen.
Dilution Range	Optimal working dilutions should be determined experimentally by the investigator Suggested starting dilutions are as follows: IF 1:2000
Formulation	Liquid in PBS pH7.4, containing 0.02% Sodium Azide and 50% Glycerol.
Isotype	IgG
Storage	Stable for one year at -15°C to -25°C from date of shipment. For maximum recovery of product, centrifuge the original vial after thawing
Instruction	and prior to removing the cap. Aliquot to avoid repeated freez

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	
Immunogen Region	
Specificity	Flag-Tag Monoclonal Antibody (2C5) AbFluor™ 405 Conjugated specially designed for your Immunofluorescence analysis.
Immunogen Sequence	