

Anti-ECFP-Tag antibody [8F6] {HRP} (STJA0006052)

STJA0006052

GENERAL INFORMATION

| | |
|--------------------------|--|
| Product Type | Primary antibodies |
| Short Description | Mouse monoclonal antibody anti-ECFP-Tag-Tag is suitable for use in Western Blot research applications. |
| Applications | WB |
| Host/Source | Mouse |
| Reactivity | Species independent |

PRODUCT PROPERTIES

| | |
|-----------------------|--|
| Clonality | Monoclonal |
| Clone ID | 8F6 |
| Concentration | 1 mg/mL |
| Conjugation | HRP |
| Purification | The antibody was purified using affinity-chromatography using specific immunogen. |
| Dilution Range | Optimal working dilutions should be determined experimentally by the investigator Suggested starting dilutions are as follows: WB 1:5000 |
| Formulation | Liquid in PBS pH7.4, containing 0.02% Sodium Azide and 50% Glycerol. |
| Isotype | IgG |
| Storage | Stable for one year at -15°C to -25°C from date of shipment. For maximum recovery of product, centrifuge the original vial after thawing |
| Instruction | and prior to removing the cap. Aliquot to avoid repeated freezing |

TARGET INFORMATION

| | |
|---------------------------|---|
| Gene ID | |
| Gene Symbol | |
| Uniprot ID | |
| Immunogen | |
| Immunogen Region | |
| Specificity | ECFP-Tag Monoclonal Antibody (8F6) HRP Conjugated, specially designed for your Western blot analysis. |
| Immunogen Sequence | |