

Anti-GAPDH antibody [2B8] {AbFluor™ 350} (STJA0006053)

STJA0006053

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Glyceraldehyde-3-phosphate dehydrogenase is suitable for use in Western Blot, Immunofluorescence and Immunohistochemistry research applications.
Applications	WB/IF/IHC
Host/Source	Mouse
Reactivity	Human/Rat/Mouse/Monkey/Dog/Chicken/Hamster/Rabbit/Pig/Sheep/Insect/Yeast

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	2B8
Concentration	1 mg/mL
Conjugation	AbFluor™ 350
Purification	The antibody was purified using affinity-chromatography using specific immunogen.
Dilution Range	Optimal working dilutions should be determined experimentally by the investigator Suggested starting dilutions are as follows: IHC 1:200 IF 1:200
Formulation	Liquid in PBS pH7.4, containing 0.02% Sodium Azide and 50% Glycerol.
Isotype	IgG1
Storage	Stable for one year at -15°C to -25°C from date of shipment. For maximum recovery of product, centrifuge the original vial after thawing
Instruction	and prior to removing the cap. Aliquot to avoid repeated freezi

TARGET INFORMATION

Gene ID	2597
Gene Symbol	GAPDH
Uniprot ID	G3P_HUMAN
Immunogen	
Immunogen Region	
Specificity	GAPDH Monoclonal Antibody (2B8) AbFluor™ 350 Conjugated specially designed for your Immunofluorescence analysis.
Immunogen Sequence	