

Anti-
HIST1H3A/HIST1H3B/HIST1H3C/HIST1H3D/HIST1H3E/HIST1H3F/HIST1H3G/HIST1H3H/HIST1H3I/HIST1H3J/HIST2H3A/HIST2H3C/HIST2H3D
antibody {AbFluor™ 350} (STJA0006091)
STJA0006091

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Histone H3.1/Histone H3.2/Histone H3.3 is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Mouse
Reactivity	Zebrafish

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	AbFluor™ 350
Purification	The antibody was purified using affinity-chromatography using specific immunogen.
Dilution Range	Optimal working dilutions should be determined experimentally by the investigator Suggested starting dilutions are as follows: IHC 1:50-300 IF 1:200-400
Formulation	Liquid in PBS pH7.4, containing 0.02% Sodium Azide and 50% Glycerol.
Isotype	IgG1
Storage Instruction	Stable for one year at -15°C to -25°C from date of shipment. For maximum recovery of product, centrifuge the original vial after thawing and prior to removing the cap. Aliquot to avoid repeated freezi

TARGET INFORMATION

Gene ID	3020/3021 8350/8351/8352/8353/8354/8355/8356/8357/8358/8968
Gene Symbol	H3-3A.H3-3B H3C1.H3C2.H
Uniprot ID	H33_HUMAN H31_HUMAN H32_HUMAN <
Immunogen	
Immunogen Region	
Specificity	Histone H3 Monoclonal Antibody AbFluor™ 350 Conjugated specially designed for your Immunofluorescence analysis.
Immunogen Sequence	