

Anti-VSV-G-Tag antibody (501-511) [8D6] {AbFluor™ 405} (STJA0006172)

STJA0006172

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-VSV-G-Tag-Tag (501-511) is suitable for use in Western Blot, Immunoprecipitation and Immunofluorescence research applications.
Applications	WB/IP/IF
Host/Source	Mouse
Reactivity	Species independent

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	8D6
Concentration	1 mg/mL
Conjugation	AbFluor™ 405
Purification	The antibody was purified using affinity-chromatography using specific immunogen.
Dilution Range	Optimal working dilutions should be determined experimentally by the investigator Suggested starting dilutions are as follows: IF 1:1000
Formulation	Liquid in PBS pH7.4, containing 0.02% Sodium Azide and 50% Glycerol.
Isotype	IgG
Storage	Stable for one year at -15°C to -25°C from date of shipment. For maximum recovery of product, centrifuge the original vial after thawing
Instruction	and prior to removing the cap. Aliquot to avoid repeated freezi

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	A synthetic peptide corresponding to aa 501-511 of vesicular stomatitis virus glycoprotein (VSV-G). This vector is a pseudo lentiviral vector.
Immunogen Region	501-511
Specificity	VSV-G-Tag Monoclonal Antibody (8D6) AbFluor™ 405 Conjugated specially designed for your Immunofluorescence analysis.
Immunogen Sequence	