

Anti-CDH3 antibody (450-550) [PTR2204] (STJA0006250)

STJA0006250

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Cadherin-3 (450-550) is suitable for use in Western Blot, Immunohistochemistry and ELISA research applications.
Applications	WB/IHC/ELISA
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	PTR2204
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was purified using affinity-chromatography using specific immunogen.
Dilution Range	WB 1:500-2000 IHC 1:50-200
Formulation	Liquid in PBS containing 50% Glycerol and 0.02% Sodium Azide.
Isotype	
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	1001
Gene Symbol	CDH3
Uniprot ID	CADH3_HUMAN
Immunogen	Synthesized peptide derived from human P-cadherin. AA range: 450-550
Immunogen Region	450-550
Specificity	This antibody detects endogenous levels of P-cadherin at Human; Predict react with mouse, rat
Immunogen Sequence	