

Anti-ATP1A1 antibody (150-250) [PTR2543] (STJA0006258)

STJA0006258

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Sodium/potassium-transporting ATPase subunit alpha-1/K (+) ATPase alpha-1 subunit (150-250) is suitable for use in Immunohistochemistry and ELISA research applications.
Applications	IHC/ELISA
Host/Source	Mouse
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	PTR2543
Concentration	
Conjugation	Unconjugated
Purification	Purified using Protein G purification.
Dilution Range	IHC-p 1:50-200
Formulation	Liquid in PBS pH7.4 containing 50% Glycerol and 0.03% Proclin 300.
Isotype	IgG3k
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	476 480 477 ATP1A1 ATP1A4 AT1A1_HUMAN AT1A4_HUMAN AT1A2_HUMAN
Immunogen	Synthesized peptide derived from human Sodium Potassium ATPase AA range: 150-250
Immunogen Region	150-250
Specificity	This antibody detects endogenous levels of Sodium Potassium ATPase at Human, Mouse, Rat
Immunogen Sequence	



Human kidney tissue was stained with Anti-Sodium Potassium ATPase (PTR2543) Antibody

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
 St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081