

Anti-APOB antibody (STJA0006288)

STJA0006288

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Camel monoclonal antibody anti-Apolipoprotein B-100 is suitable for use in Flow Cytometry and ELISA research applications.
Applications	FC/ELISA
Host/Source	Camel
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified nickel column using his-tag in antibody
Dilution Range	ELISA 1:10000-20000
Formulation	Liquid in PBS
Isotype	
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	338
Gene Symbol	APOB
Uniprot ID	APOB_HUMAN
Immunogen	Recombinant protein of full length human ApoB
Immunogen Region	
Specificity	This antibody detects endogenous levels of ApoB
Immunogen Sequence	