

**Anti-LCN2 antibody (STJA0006289)**

STJA0006289

**GENERAL INFORMATION**

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Camel monoclonal antibody anti-Neutrophil gelatinase-associated lipocalin is suitable for use in Flow Cytometry and ELISA research applications.
<b>Applications</b>	FC/ELISA
<b>Host/Source</b>	Camel
<b>Reactivity</b>	Human/Mouse

**PRODUCT PROPERTIES**

<b>Clonality</b>	Monoclonal
<b>Clone ID</b>	
<b>Concentration</b>	1 mg/mL
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	The antibody was affinity-purified nickel column using his-tag in antibody
<b>Dilution Range</b>	ELISA 1:10000-20000
<b>Formulation</b>	Liquid in PBS
<b>Isotype</b>	
<b>Storage Instruction</b>	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

**TARGET INFORMATION**

<b>Gene ID</b>	3934
<b>Gene Symbol</b>	LCN2
<b>Uniprot ID</b>	NGAL_HUMAN
<b>Immunogen</b>	Recombinant protein of full length human NGAL.
<b>Immunogen Region</b>	
<b>Specificity</b>	This antibody detects endogenous levels of NGAL.
<b>Immunogen Sequence</b>	